

C2004 400ml Versatile centrifuge

"entry level...excellence"

- Speed** 500-4000Rpm
in 10 Rpm steps
- Max Rcf** 2200G
- Timer** 0-99Minutes
and constant
- Dims HWD** 27x33x46Cms
- Weight** 16.5Kilos
without rotor



Rotors Available



BR5431

max speed 4000 Rpm 2,200G

24 x 15ml tube, fixed angle rotor
Tube size maximum H120 x W17 mm
5ml and 10ml

Reducers available



BR5434

max speed 4000 Rpm 2,200G

8 x 50ml tube, fixed angle rotor
Tube size maximum H115 x W30mm
15ml and 25ml

Reducers available



BR5505

max speed 4000 Rpm 2,200G

4 x 100ml max swing out rotor with rings
Requires tube adaptors as below
Sealed lids available B5419



TUBE ADAPTORS for BR5505

Range from 0.5ml to 100ml max

See Page 14/15 for full details