

C2015 Micro high speed centrifuge

“speedy little mover”

Speed 500-15000Rpm
in 10 Rpm steps

Max Rcf 22000G

Timer 0-99Minutes
and constant

Dims HWD 23x23x35Cms

Weight 12Kilos
without rotor



Rotors Available



BR5424

max speed 15000 Rpm 22,000G
24 x 1.5-2.2ml tube, fixed angle rotor
Tube size maximum H40 x W11mm
Cap size dependant
0.2ml - 0.4ml and 0.5ml

Reducers available



BR5436

max speed 15000 Rpm 22,000G
36 x 0.5ml tube, fixed angle rotor
Tube size maximum H40 x W8mm
Cap size dependant



BR5448

max speed 15000 Rpm 22,000G
48 x 0.2ml tube, fixed angle rotor
Tube size maximum H40 x W6mm
Cap size dependant



BR5494

max speed 15000 Rpm 22,000G
4 x 8 place PCR strip (0.2ml)