

C2041 1 Litre Centrifuge

“from the family of workhorses”

- Speed 500-4000Rpm in 10 Rpm steps
- Max Rcf 2600G
- Timer 0-99Minutes and constant
- Dims HWD 28x40x50Cms
- Weight 22Kilos without rotor



Rotors Available



BR5431

max speed 4000 Rpm 2,200G
24 x 15ml tube, fixed angle rotor
Tube size maximum H120 x W17 mm
5ml and 10ml

Reducers available



BR5434

max speed 4000 Rpm 2,200G
8 x 50ml tube, fixed angle rotor
Tube size maximum H115 x W30mm
15ml and 25ml

Reducers available



BR5510

max speed 4000 Rpm 2,600G
Swing out rotor 4 x 250ml max
Requires tube adaptors as below

B5250

Buckets required set of 4 (shown with Rotor)
Sealed lids available B5419



TUBE ADAPTORS for BR5510

Range from 0.5ml to 250ml max

See Page 14/15 for full details