

C2043 3 Litre work horse

"Big is best"

Speed 500-4000Rpm
in 10 Rpm steps

Max Rcf 3600G

Timer 0-99Minutes
and constant

Dims HWD 40x60x63Cms

Weight 43Kilos
without rotor



Rotors Available



BRK5553

max speed 4000 Rpm 3,600G
4 x 750ml max swing out rotor
Requires B3000 buckets as below



B3000

max speed 4000 Rpm 3,600G
Buckets for above set of 4
Sealed lids available B5319

ADAPTORS

for B3000
See Page 14/15 for full details

Oil Adaptors available - Please call.