

CR2000

Multi purpose & versatile rotor recognised centrifuge

"Top selling all rounder"

Speed 500-12000Rpm
in 10 Rpm steps

Max Rcf 13500G

Timer 0-99Minutes
and constant

Dims HWD 27x33x46Cms

Weight 16.5Kilos
without rotor



Rotors Available



BRK5505

max speed 4000 Rpm 2,200G

4 x 100ml max swing out rotor with rings
Requires tube adaptors as below
Sealed lids available B5419



TUBE ADAPTORS for BRK5505

Range from 0.5ml to 100ml max

See Page 14/15 for full details



BRK5434

max speed 4000 Rpm 2,200G

8 x 50ml tube, fixed angle rotor
Tube size maximum H115 x W30mm
15ml and 25ml

Reducers available



BRK5431

max speed 4000 Rpm 2,200G

24 x 15ml tube, fixed angle rotor
Tube size maximum H120 x W17 mm
5ml and 10ml

Reducers available

Rotors - fixed angle



BRK5401

max speed 12000 Rpm 13,500G
24 x Capillary tubes
Tube size maximum H75 x W1.6mm
Also 12 x 1.5-2.0ml tube
Tube size maximum H40 x W11mm
Cap size dependant.



BRK5424

max speed 12000 Rpm 13,500G
24 x 1.5-2.2ml tube, fixed angle rotor
Tube size maximum H40 x W11mm
Cap size dependant.

Reducers available

0.2ml - 0.4ml and 0.5ml



BRK5436

max speed 12000 Rpm 13,500G
36 x 0.5ml tube, fixed angle rotor
Tube size maximum H40 x W8mm
Cap size dependant

BRK5448

48 x 0.2ml tube, fixed angle rotor
Tube size maximum H40 x W6mm.



BRK5494

max speed 12000 Rpm 13,500G
4 x 8 place strip PCR rotor (0.2ml)



BRK5412

max speed 6000 Rpm 3,800G
12 x 15ml tube, fixed angle rotor
Tube size maximum H120 x W17 mm

Reducers available

5ml and 10ml



BRK5406

max speed 6000 Rpm 3,800G
6 x 50ml Tube, fixed angle rotor
Tube size maximum H115 x W30mm

Reducers available

15ml and 25ml