

# C2006 Fixed angle up to 6000RPM

*“functional simplicity”*

- Speed** 500-6000Rpm  
in 10 Rpm steps
- Max Rcf** 3800G
- Timer** 0-99Minutes  
and constant
- Dims HWD** 27x33x46Cms
- Weight** 16.5Kilos  
without rotor



## Rotors Available



**BR5412**

max speed 6000 Rpm 3,800G  
12 x 15ml tube, fixed angle rotor  
Tube size maximum H120 x W17 mm



**RM10**

ADAPTORS for BR5412 (above) 5-10ml  
Tube size maximum H100 x W13 mm



**BR5406**

max speed 6000 Rpm 3,800G  
6 x 50ml Tube, fixed angle rotor  
Tube size maximum H115 x W30mm



**RM25**

ADAPTORS for BR5406 (above) 25ml  
Tube size maximum H100 x W25 mm

**RM15**

ADAPTORS for BR5406 (above) 15ml  
Tube size maximum H120 x W17 mm