

K2045 *"Designed with your every need in mind"*

Speed 500-15000Rpm
in 1 Rpm steps

Max Rcf 22,000G

Timer 0-999Minutes
+ hold in
1 second steps

Dims HWD 27x33x46Cms

Weight 17Kilos
without rotor



Rotors - swing out



BRK5505

max speed 4000 Rpm 2,200G
4 x 100ml max swing out rotor
with rings
Requires tube adaptors as below
Sealed lids available B5419



TUBE ADAPTORS for BR5505

Range from 0.5ml to 100ml max

See Page 14/15 for full details

Rotors - fixed angle



BRK5324

max speed 6000 Rpm 4,800G
24 x 15ml tube, fixed angle rotor
Tube size maximum H120 x W17 mm
Reducers available 5ml and 10ml



BRK5308

max speed 6000 Rpm 4,800G
8 x 50ml tube, fixed angle rotor
Tube size maximum H115 x W30mm
Reducers available 15ml and 25ml



BRK5494

Rotors - fixed angle

max speed 15000 Rpm 22,000G
4 x 8 place strip PCR rotor (0.2ml)



BRK5424

max speed 15000 Rpm 22,000G
24 x 1.5-2.2ml tube, fixed angle rotor
Tube size maximum H40 x W11mm
0.2ml - 0.4ml and 0.5ml

Reducers available



BRK5436

max speed 15000 Rpm 22,000G
36 x 0.5ml tube, fixed angle rotor
Tube size maximum H40 x W8mm



BRK5448

max speed 15000 Rpm 22,000G
48 x 0.2ml tube, fixed angle rotor
Tube size maximum H40 x W6mm