










Adaptors for Swing out rotors

Adaptors

Adaptors

Adaptors

Tube type	Shape	Capacity	Size	To fit rotor BR5505		To fit rotor BR5510		To fit rotor BR5553		
				Part nos	Nos of tubes Per rotor	Need buckets Part nos	Nos of tubes Per rotor	Need buckets Part nos	Nos of tubes Per rotor	
	Micro with cap	point	0.5ml	8x20	AM505	40	AM605	40	AM805	120
	Micro with cap	point	1.5ml	11x38	AM520	28	AM620	28	AM820	84
	Micro with cap	point	2.0ml	11x38	AM520	28	AM620	28	AM820	84
	Micro with cap	point	0.2ml	6x20	AM502	48	AM602	64	AM802	168
	Micro with cap	point	0.4ml	6x30	AM504	48	AM604	64	AM804	168
	Plain No cap	round	1ml	6x45	AR501	48	AR601	48	AR801	144
	Plain No cap	round	3ml	10x60	AR503	32	AR603	32	AR803	96
	Plain No cap	round	5ml	12x75	AR505	32	AR605	32	AR805	96
	Plain No cap	round	6ml	12x82	AR505	32	AR605	32	AR805	96
	Plain No cap	round	7ml	12x100	AR507	32	AR607	32	AR807	96
	Plain No cap	round	9/10ml	14x100	AR509	28	AR609	28	AR809	96
	Plain No cap	round	15ml	17x100	AR515	28	AR615	28	AR815	52
	Plain No cap	round	25ml	24x100	AR525	12	AR625	12	AR825	28
	Plain No cap	round	50ml	34x100	AR550	4	AR650	8	AR850	16
	Plain No cap	round	100ml	45x100	AR5100	4	AR6100	4	AR8100	8
	Plain No cap	round	150ml	52x100	N/A	N/A	AR6150	4	AR8150	4
	Plain No cap	round	250ml	62x100	N/A	N/A	BUCKET	4	AR8250	4
	Falcon with cap	Point	15ml	17x120	AF515	8	AF615	16	AF815	48
	Falcon with cap	Point	50ml	29x115	AF550	4	AF650	4	AF850	20
	Falcon with cap	Point	175ml	61x118	N/A	N/A	AF6175	4	AF8175	4
	Corning with cap	Point	250ml	60x172	N/A	N/A	N/A	N/A	AF8250	4
	Corning with cap	Point	500ml	98x148	N/A	N/A	N/A	N/A	AF8500	4
	Falcon with cap	Square	12ml	17x100	AFS512	8	AFS612	16	AFS812	48
	Falcon with cap	Square	25ml	25x90	AFS525	4	AFS625	8	AFS825	28
	Falcon with cap	Square	30ml	25x110	AFS530	4	AFS630	8	AFS830	28
	Falcon with cap	Square	50ml	29x115	AFS550	4	AFS650	4	AFS850	20
	Nalgene/Oakridge	Round	10ml	16x80	ANO510	8	ANO610	16	ANO810	64
	Nalgene/Oakridge	Round	30ml	26x95	ANO530	4	ANO630	8	ANO830	28
	Nalgene/Oakridge	Round	50ml	29x107	ANO550	4	ANO650	8	ANO850	24
	Nalgene/Oakridge	Round	85ml	38x106	ANO585	4	ANO685	4	ANO885	12
	Nalgene/Oakridge	Flat	250ml	62x130	N/A	N/A	ANO6250	Bucket	ANO8250	4
	Nalgene/Oakridge	Flat	750ml	98x153	N/A	N/A	N/A	N/A	ANO8750	4
	Monovette	Square	1.1-1.4ml	8x82	AM5014	48	AM6014	48	AM8014	96
	Monovette	Square	2.7-3ml	11x82	AM503	48	AM603	48	AM803	76
	Monovette	Square	2.6-2.9ml	13x81	AM529	40	AM629	40	AM829	76
	Monovette	Square	4.5-5ml	11x108	AM503	48	AM603	48	AM803	76
	Monovette	Square	7.5-8.2ml	13x106	AM579	40	AM679	40	AM879	76
	Monovette	Square	4-5.5ml	15x92	AM550	28	AM650	28	AM850	64
	Monovette	Square	7.5-8.2ml	15x108	AM550	28	AM650	28	AM850	64
	Monovette	Square	9-10.0ml	16x108	AM590	28	AM690	28	AM890	64
	Vacutainer	Round	1.6-5.0ml	13x75	AV516	40	AV616	40	AV816	76
	Vacutainer	Round	4-7.0ml	13x100	AV550	40	AV650	40	AV850	76
	Vacutainer	Round	4-7.0ml	16x75	AV570	28	AV670	28	AV870	64
	Vacutainer	Round	8.0ml	16x125	AV580	28	AV680	28	AV880	64
	Vacutainer	Round	8.5-10.0ml	16x100	AV580	28	AV680	28	AV880	64